

MEDILINK^{esi} 2010+ (R1.00)
Feature List

Patrick Sims
Developer - Medilink
2009/12/07

Update 2010/07/08

Table Of Contents

Table Of Contents.....	ii
Introduction.....	1
HICAPS Integration.....	2
ICTK Integration.....	3
Medilink Scanner/Document Manager.....	4
Other Features.....	5
Next Available Appointment For Practice.....	5
Next Appointment Label.....	5
2Hippo SMS Reminders.....	6
Convert Quote To Invoice.....	6
Reverse Receipts.....	6
Different Versions.....	6
Healthcare Identifiers.....	6
Prompt On Close.....	7
Update Provider Leave.....	7
Invoice Dialog.....	7
Revamp.....	7
Edit Item.....	8
Patient Appointment Log.....	8
Reports.....	10
Provider Holidays Report.....	10

Introduction

This document was written by Patrick Sims to detail the various new features in the Medilink^{ESI} 2010+ R1 release.

HICAPS Integration

Medilink has gone through an integration procedure with HICAPS. These key features of the integration are:

- Medicare Easyclaim
- Integrated EFTPOS
- HICAPS Claiming


What this means for our clients is that they can now perform these tasks through the Medilink interface with the benefits of:

- No double keying in to EFTPOS and software.
- Much quicker than keying data in to an EFTPOS terminal.
- Reports and other information all available through Medilink.

This integration requires that you obtain a HICAPS/NAB device.

Please see the Medilink HICAPS documentation and/or contact Medilink for more information.

A full user guide is available for users and installers.



ICTK Integration

Medilink has also integrated with iCTK. This product is a joint project between iSOFT and Touch Networks. It provides the following modules:

- Medicare Easyclaim
- Integrated EFTPOS


The benefits of using this integration system are:

- No double keying in to EFTPOS and software.
- Much quicker than keying data in to an EFTPOS terminal.
- Reports and other information all available through Medilink and iCTK software.

ICTK is unique in that it will support a variety of bank terminals.

Please see the Medilink iCTK documentation and/or contact Medilink for more information.


A full user guide is available for users and installers.



Medilink Scanner/Document Manager

Medilink has produced a new scanning/document manager module. This allows for scanning through any TWAIN compatible scanner, (most scanners will come with TWAIN drivers).


Additionally, clients on existing Redmap or other scanning modules can import their existing data in to the new system.



Other Features

Next Available Appointment For Practice

We have extended the Next Available Appointment feature to allow selection of a practice:



Next Appointment Label

There are now two additional options when printing out next appointment labels, found in 'General Options'.



You can elect to turn on/off the 'headers' in the label, (i.e., 'Name:', 'Seeing:', etc.). You can also elect to print your registered address in the label, (please note, you should turn this off if you print labels for other practices, as the address will be incorrect).

Your Next Appointment:

Name:

Seeing:


Location:

Phone Number:

11:30am
 Wednesday July 14, 2010
 Mr Leslie EELES
 Dr Miriam Gadney
 This Practice 31 Street St
 00 0000 0000

2Hippo SMS Reminders

The existing SMS reminders system has been revamped and the existing bugs ironed out. This still carries messages through 2Hippo, and allows for 2 way patient confirmations.



Please contact Medilink for more information.

Convert Quote To Invoice

In the 'Copies' dialog you can now convert quotes to an invoice. This is handy when combined with the new 'edit item' feature.

Reverse Receipts

Reverse receipts now detail the type of EFTPOS transaction, (SAV/CHQ/CR).

Different Versions

You will now not be given the choice to run Medilink if it is detected that a different version is running on the server.


Healthcare Identifiers

New fields have been added for the NEHTA standard healthcare identifiers, (National E-Health Transition Authority, a government agency). These are universal identifiers for patients, providers and practices, (organisations), that provide everyone with standard methods of identification. This in turn will improve the inoperativeness of medical systems, and have a flow on effect for the speed and quality of patient care.

These will be unlocked later in the year when added functionality, (including using web services to obtain patient/provider/referrer healthcare identifiers), and the standard is formally adopted, (they are currently the combination of two existing international standards, ISO 7812-1:2006 and HL7 OIDs).

Prompt On Close

There is now an option that will confirm that you want to quit.



Update Provider Leave

Medilink can now update provider leave dates in the appointment book interface at regular intervals.

Invoice Dialog

Revamp

The standard invoice dialog has been revamped, (to enable the other types of billing that it facilitates, such as Easyclaim).

Private Invoicing (HIC Online)

Account: Mr Leslie EELES Claimant: EELES Leslie Mr 2298 43130 11
 Invoiced Patient: Mr Leslie EELES 2298 43130 11 Payee Provider: Dr. M. Gadney 2107501Y
 Searched Patient: Mr Leslie EELES 2298 43130 11 Referring Provider: Dr Brenda Reed 2054781W
 Provider: Dr Miriam Gadney Referral Issue Date: 22/04/2010 Referral Period: Standard 12 months
 Practice: This Practice 2107501Y
 Location: This Practice
 Practice Type: S - Specialist

Service Date	Time (24hr)	Item Number	Description	Send Service Text	GST Included	Fee	Code
12/07/2010							

Buttons: Calendar, Dur (min), List, History, Multiple, Hide, Edit, Delete, Add GST

Service	Item	Description	Code	Fee	GST

Print Options: Invoice Claim Form (HIC Online) Use Crystal Extras CutePDF Writer Total \$
 Use Med Addr 3rd P Bulk Billing Rebate STS Amt (Hosp) Deferred Billing

HIC Online
 Claims will be sent interactively
 Setup: Links | HIC Online | Configure

Buttons: Prepay, Quote, Cancel, Help


It now shows more information about the invoice/claim, GST, and the code for example, is now in line with the items that you add.

Edit Item

You can now edit items that have previously been added. Highlighting an item that you have added, then pressing the button will remove that item from the list, and put it back in the fields up above. You then make the changes that you want, and press 'Add' again. This comes in handy for converting quotes to invoices.

Patient Appointment Log

The patient appointment log has been revamped to allow for searching by appointment date or user action, and to filter out the status changes and only show add new/copy/move/delete appointments.



The image shows a 'User Report Search' dialog box. At the top, there is a title bar with the text 'User Report Search' and a close button. Below the title bar, there is a text input field containing 'User Activity Report - Appointments' and a checkbox labeled 'Use In Search'. The main area of the dialog contains two radio buttons: 'Date of User Action' (unselected) and 'Appointment Date' (selected). Below these are two rows for date selection: 'Start Date' and 'End Date', both with text boxes containing '14/07/2010' and 'Calendar' buttons. To the right of each date row is a checked checkbox. Below the dates are two text input fields: 'User Name:' (empty) and 'Patient Name:' (containing 'SMITH'). To the right of each name field is a checkbox, with the 'Patient Name' checkbox checked. At the bottom of the main area is a checkbox labeled 'Show Add New/Copy/Move/Delete Appointments Only?' which is checked. At the bottom of the dialog are two buttons: 'OK' with a green checkmark icon and 'Cancel' with a red X icon.

Reports

Provider Holidays Report

There is now a report to list provider leave dates, (holidays):

Medilink System Report
Provider Holidays - All Holidays



Report Date 8/7/10 2:02:37 PM

<u>Provider</u>	<u>Date From</u>		<u>Start Time</u>	<u>Date To</u>		<u>Finish Time</u>	<u>Days</u>
<u>Gadnev, Miriam Dr</u>	2010/01/01	Friday	800	2010/01/05	Tuesday	2000	5